



Safe Pick® Sharps Pouch

Code: HA1887FLP



For your own safety read these instructions fully

All Safe Pick Pouch items can be ordered as replacements—please request our Safe Pick Spares Order Form from Head Office.

	Qty	Description
Pair of Disposable Nitrile Gloves Code: HX01890	1	These Nitrile gloves surpass water retention and viral penetration tests, hence they help protect against penetration of viruses to the skin. They do not trigger Latex allergies, are comfortable to wear and easy to remove without self-contamination, following the removal procedure below.
Pair of Stainless Steel Forceps Code: HX01920	1	Grooves at the forcep ends help to pick up a needle (even with a syringe attached) and a ratchet between the handles locks the forceps round an object until the time of releasing the needle into the isolation cylinder. The steel forceps give more grip to extract needles. Take care not to release accidentally as this may cause the needle to fly out. The forceps are ideal for instances where it is not practical to use the Safe Pick Litterpicker.
Safe Pick Litterpicker Code: HX01880	1	The essential in safety equipment, the folding Safe Pick Litterpicker allows a safe working distance between the user and needle/syringe. The user must ensure the needle/syringe is secured firmly between Litterpicker jaws whilst it is being transferred to the Needle Isolation Cylinder.
Needle Isolation Cylinder 0.1ltr Code: HX01900	1	Wherever practical, lie the Cylinder on the ground adjacent to the needle/syringe. DO NOT HOLD THE CYLINDER WHILST TRANSFERING THE NEEDLE INTO IT. The lid can be pushed close to prevent the needle falling out, and if you push the lid hard closed it will lock. Only do this when you are ready to dispose of the cylinder and contents; this must be done via a Licensed Clinical Waste disposal company. Always treat an occupied, a full or partially full cylinder as Clinical Waste. NB a Sharps Disposal Container 0.6ltr is also available on request.
Disinfectant Spray Code: HX01905	1	Essential for the thorough cleaning and disinfection of hands and equipment (e.g. Litterpicker) and the contaminated area.
Disinfectant Wipes Code: HX01910	5	Essential for the thorough cleaning and disinfection of hands and equipment (e.g. Litterpicker) and the contaminated area.
Bio Hazard Bag Code: N/A	4	For the safe containment of items for disposal.

Instructions For Use

1. First, put on a pair of the disposable Nitrile Gloves.
2. Have the Isolation Cylinder open and ready (laying on the ground adjacent to the needle - if practical).
3. Use the Safe Pick Litterpicker or forceps to remove the needle from the area and release into the isolation cylinder.
4. Tightly close the lid and dispose of in a yellow clinical waste sack.
5. Thoroughly clean and disinfect any re-usable equipment.
6. Seal the yellow clinical waste sack and dispose of (or even transported within) an appropriately sized sharps bin.
7. Carefully remove the disposable gloves and also put them into the yellow clinical waste bag.
8. Discard all waste into the yellow clinical waste sack with all packaging.
9. Use the hand wipes to thoroughly cleanse your hands.
10. Seal the yellow clinical waste sack and dispose of (or even transported within) an appropriately sized sharps bin.
11. As soon as possible, wash your hands in hot, soapy water.

Important Notes

Clinical waste must NEVER be disposed of in normal domestic or trade refuse. These instructions are suggestions for use. Please ensure you comply with your company Health & Safety regulations at all times.

Contact Details:

Postal: Bromyard Road, Ledbury, Herefordshire, HR8 1NS

Tel: 01531 635678; Fax: 01531 635670

Email: sales@helpinghand.co.uk - www.litterpickersdirect.com



#1 for litter clearance hand tools